



### General Medical Bed Head Unit



We are the manufacturer of Zenith factory, which specializes in selling cheap and durable products made in China. Our DFTT® General Medical Bed Head Unit adopts top-notch materials to meet the basic needs of wards and is exported to many countries. It is convenient and fast to install and has no harsh requirements for users. At the same time, General Medical Bed Head Unit can be equipped with multiple medical tools. We would like to be your partner in China.

## Product Description

DFTT® **General Medical Bed Head Unit** is strong and easy to use. **General Medical Bed Head Unit** is made of A6061 material and has a irregular shape. We have a wide variety of Medical Bed Head Unit, this product can meet the basic needs of ordinary ward, simple without losing its basic function, can also be configured separately related equipment, can be customized, to meet all the needs of customers.

### Product Parameter

Product name	size	material	Gas type	color
<b>General Medical Bed Head Unit</b>	118*30*100CM	A6061	All	white customized

### General Medical Bed Head Unit Feature And Application





Medical Bed Head Unit is often used in the bedside of the ward, can be connected to a variety of equipment, to ensure patient safety. Medical Bed Head Unit is characterized by safety and reliability, which is conducive to doctors' timely treatment of patients. The overall design of the Medical Bed Head Unit is reasonable, saving space to the maximum extent, and providing a comfortable environment for patients. It integrates practicability, security and other performance, and is easy to maintain. In the aspect of after-sale, we will solve and correct all quality problems. Medical Bed Head Unit Details.



[www.clear-medical.com](http://www.clear-medical.com)



Ningbo Qingjie Medical Devices Co.,Ltd  
Phone:+86-13924646868  
E-mail: dftchina@163.com



DFTT® Medical Bed Head Unit is composed of socket, switch and terminal jack. Medical Bed Head Unit is easy to install and firm. Medical Bed Head Unit is sealed and connected between plates. The biggest advantage of the Medical Bed Head Unit is that it simplifies the traditional bedside shelf and integrates modern technology. In view of the product instructions and precautions, there is no special matters, is to handle with care.