

Japanese Standard Gas Outlets



We are a medical device company in China, and we are also a manufacturer, we can produce the products you need by ourselves. DFTT® Japanese Standard Gas Outlets is one of our most popular products, which has been exported to many foreign countries. Japanese Standard Gas Outlets are affordable, safe and reliable, high quality, which is your best choice. We are looking forward to becoming your partner in China.

Product Description

DFTT® **Japanese Standard Gas Outlets** are made of copper material, with strong and stable structure and excellent quality, and can also be customized according to individual needs of customers. **Japanese Standard Gas Outlets** are specialized in design, excellent in workmanship and complete in function. They can match and exchange with the original ones and realize the upgrading and upgrading of new Outlets. They are excellent

Product Parameter

Product name	Shape	material	size	color
Japanese Standard Gas Outlets	irregular shape	Copper	7*6*6CM	Golden

Japanese Standard Gas Outlets Feature And Application





Japanese Standard Gas Outlets have passed the quality system certification. They are safe and reliable, and can be distributed and exchanged with original products to realize upgrading and upgrading. **Japanese Standard Gas Outlets** are compact and convenient for connecting many kinds of Gas pipelines and equipped to form a whole bed unit. In terms of layout, **Japanese Standard Gas Outlets** complement each other with the whole cooperation, and the color and style are in line with the design of the ward, giving the patients a warm feeling.

Japanese Standard Gas Outlets Details



Japanese Standard Gas Outlets built-in access valve, can avoid disassembly panel with air maintenance, easy installation and maintenance, not easy to cause damage. **Japanese Standard Gas Outlets** distinguish Gas types by using different pin holes in different locations and have a single configuration to prevent different Gas Outlets from being misconnected and reduce the occurrence of safety accidents. **Japanese Standard Gas Outlets** integrated all-metal material, durable, and passed the factory air tightness test. The surface of the Gas Outlets has been equipped with advanced processing technology, which has passed the quality system Standard and 20000 plug and pull fatigue tests.